


KOMET® KONNECT



T-Post™ System Reconstruction made simple

With posts fabricated of pure titanium, the KOMET® T-Post™ system simplifies reconstruction procedures and enhances reliable, restorative outcome. Posts are available in a variety of head designs and are offered in several sizes and diameters and in four shank lengths (6, 7, 9, and 13 mm). The T-Post™ cylindrical-head post permits reconstruction of partially destroyed or multi-rooted teeth or in cases requiring the placement of multiple posts. The T-Post™ X post has a “step-down” design on its coronal portion to simplify the reconstruction of severely damaged teeth, and the T-Post™ X ELO post is specially designed for the reconstruction of severely destroyed teeth with long clinical crowns. Both the T-Post™ X and T-Post™ X ELO posts have individually segmented heads to provide reliable, secure retention of restorative materials. Each post in the KOMET® T-Post™ system is available in active (threaded) or passive versions to correspond to different root requirements.




Kit 4122

Retail \$220.20

Buy One T-Post™ X Kit 4122 at retail
Get One T-Post™ X 10-pack* **FREE**

*T63L6, T63L7, or T63L9 • \$95 Value



“The KOMET® T-Post™ system has been the best for ease of placement, durability, and choice of sizes of any system we have used. Immediate placement with build-up and prepping for initial placements or salvaging existing bridgework has been especially beneficial.”

Terry Bawcom, DDS, Abilene, TX

“Ease of use. Reliability. Core material never breaks. Outstanding assortment of sizes and shapes. All these aspects make the T-Post™ system the only post system our group uses. I have no interest in using any other system.”

Philip V. Mazzola, DDS, Nesconset, NY

