



Q-Finishers[®]

For the highest quality finish

- Does the work of two instruments
- Saves time
- Saves money
- Special cross cut tothing for working on composite and natural tooth

A world-wide leader
in rotary instruments

888-566-3887
www.komet-usa.com



Product characteristics and advantages

With our special cross cut toothing designed to work on composites and natural tooth, the Q-Finisher® is a real go-getter. Q-Finishers® offer an inexpensive way of trimming and finishing in fewer steps. Working on composite fillings usually required three finishing steps with normal, fine and extra-fine finishing instruments. With the development of an innovative toothing, the procedure can now be reduced to just two steps, using the Q-Finisher® for the first step and an Ultra-fine finishing instrument for the second step. This results in less chair time and a better result.

Save time and money with the best finisher in the business.



Advantages:

- Time saving because one step can be omitted
- Cost saving because one instrument can be omitted
- After the first finishing step, a better surface quality is achieved than previously after the second step. This is due to the cross-cut toothing specially designed for working on fillings
- The tapered figures H134Q and H135Q assure gentle finishing without damage to the gingiva, thanks to their smooth non-cutting tip

