Ordering guide | Single-Patient Use Instruments




single
patient use


Practice made perfect.




Shank types




204 RA
right angle — ∅ 2,35 mm




205 RASL
right angle · surgical length — ∅ 2,35 mm



313 FGSS
friction grip · short shank — ∅ 1,60 mm



314 FG
friction grip — ∅ 1,60 mm



316 FGSL
friction grip · surgical length — ∅ 1,60 mm

Diamond grit sizes

fine	●	46 µm
medium	●	107 µm *
coarse	●	151 µm *
super-coarse	●	181 µm *

* With some instruments the grit size may deviate from the specified value, depending on their shape and size.

Maximum permissible speed

■ = $n_{max.}$ 100 min⁻¹/*rpm*
□ = $n_{max.}$ 500 min⁻¹/*rpm*
◧ = $n_{max.}$ 800 min⁻¹/*rpm*
● = $n_{max.}$ 1.000 min⁻¹/*rpm*
○ = $n_{max.}$ 5.000 min⁻¹/*rpm*
◓ = $n_{max.}$ 6.000 min⁻¹/*rpm*
◒ = $n_{max.}$ 7.000 min⁻¹/*rpm*
◐ = $n_{max.}$ 8.000 min⁻¹/*rpm*
● = $n_{max.}$ 10.000 min⁻¹/*rpm*
○ = $n_{max.}$ 15.000 min⁻¹/*rpm*

◓ = $n_{max.}$ 16.000 min⁻¹/*rpm*
▲ = $n_{max.}$ 20.000 min⁻¹/*rpm*
△ = $n_{max.}$ 25.000 min⁻¹/*rpm*
◆ = $n_{max.}$ 30.000 min⁻¹/*rpm*
◔ = $n_{max.}$ 35.000 min⁻¹/*rpm*
◔ = $n_{max.}$ 40.000 min⁻¹/*rpm*
◇ = $n_{max.}$ 50.000 min⁻¹/*rpm*
◓ = $n_{max.}$ 60.000 min⁻¹/*rpm*
◒ = $n_{max.}$ 70.000 min⁻¹/*rpm*
◐ = $n_{max.}$ 80.000 min⁻¹/*rpm*

◑ = $n_{max.}$ 90.000 min⁻¹/*rpm*
▮ = $n_{max.}$ 100.000 min⁻¹/*rpm*
▮ = $n_{max.}$ 120.000 min⁻¹/*rpm*
▯ = $n_{max.}$ 140.000 min⁻¹/*rpm*
▯ = $n_{max.}$ 150.000 min⁻¹/*rpm*
▮ = $n_{max.}$ 160.000 min⁻¹/*rpm*
✦ = $n_{max.}$ 300.000 min⁻¹/*rpm*
$n_{max.}$ 450.000 min⁻¹/*rpm*
if not indicated

Some of the products and designations mentioned in the text are trademarked, patented or copyrighted.
The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products subject to alterations. Printing errors excepted.

As at January 2020

✆ 888-KOMET-USA · 800-223-7485

Single-Patient Use Instruments

Komet is the highest quality provider of burs for dental practices focused on individualized care that value the peace of mind that comes with having the sharpest bur, every time.


The new single-patient use line is the same quality Komet is known for - with the added benefits of confidence in always sterile tools.

Advantages:
- No sterilization required
- Ease of mind regarding patient safety
- Never dull, brand new bur is always sharp

Available in rolls of 25 for diamond and carbide burs, and rolls of 10 for all polishers, the packaging stores and dispenses the rolls while featuring a remaining inventory count on each individual package. QR codes on packaging allow for easy reordering and inventory maintenance.

Sterile, single-patient use Komet burs are quick and efficient, and especially convenient used in combination with Komet's large portfolio of multi-use rotary instruments.





single
patient use

Practice made perfect.

Tungsten Carbide Operative Carbides



¼ **H1SU.**FG.005

½ **H1SU.**FG.006

2 **H1SU.**FG.010

4 **H1SU.**FG.014

6 **H1SU.**FG.018

8 **H1SU.**FG.023 +

2 **H1SU.**FGSL.010 ∎

4 **H1SU.**FGSL.014 ∎

6 **H1SU.**FGSL.018 ∎

8 **H1SU.**FGSL.023 ∎

½ **H1SU.**RA.006 ∎

2 **H1SU.**RA.010 ∎

4 **H1SU.**RA.014 ∎

6 **H1SU.**RA.018 ∎

8 **H1SU.**RA.023 ∎

2 **H1SU.**RASL.010 ∎

4 **H1SU.**RASL.014 ∎

6 **H1SU.**RASL.018 ∎

8 **H1SU.**RASL.023 ∎



33 ½ **H2SU.**FG.006

34 **H2SU.**FG.008

35 **H2SU.**FG.010



330 **H7SU.**FGSS.008

330 **H7SU.**FG.008

331 **H7SU.**FG.010

332 **H7SU.**FG.012



330MWV **H24SU.**FG.008



245 **H245SU.**FG.008



56 **H21SU.**FG.009

57 **H21SU.**FG.010

58 **H21SU.**FG.012



57 **H21RSU.**FG.010



557 **H31SU.**FG.010

558 **H31SU.**FG.012

557 **H31SU.**FGSL.010 +

558 **H31SU.**FGSL.012 +



557S **H31SSU.**FG.010



1557S **H31RSSU.**FG.010

1558S **H31RSSU.**FG.012



169 **H23SU.**FG.009

170 **H23SU.**FG.010



169L **H23LSU.**FG.009

170L **H23LSU.**FG.010



1171 **H23RSU.**FG.012



699 **H33SU.**FG.009

700 **H33SU.**FG.010

701 **H33SU.**FG.012

702 **H33SU.**FG.016

701 **H33SU.**FGSL.012 +



701L **H33LSU.**FGSL.012 ∎

Tungsten Carbide Crown cutter



H40SU.FG.012

Tungsten Carbide Finishing Burs



H132SU.FG.008 +



H133SU.FG.010



H134SU.FG.014



H135SU.FG.014 +



7801 **H247SU.**FG.009

7802 **H247SU.**FG.010

7803 **H247SU.**FG.012



7004 **H41SU.**FG.014

7006 **H41SU.**FG.018

7008 **H41SU.**FG.023 +



7664 **H375RSU.**FG.014 +

7675 **H375RSU.**FG.016

7686 **H375RSU.**FG.018

Diamonds



H283SU.FG.012 +

H283SU.RA.012 ▮



H283KSU.FG.016



7901
H246SU.FG.009

7902
H246SU.FG.010

7903
H246SU.FG.012



H48LSU.FG.010 +

H48LSU.FG.012 +



H274SU.FG.016

H274SU.FG.018



7404
H379SU.FG.014

7406
H379SU.FG.018

7408
H379SU.FG.023

7408
H379SU.RA.023 ▮



8801SU.FG.016

8801SU.FG.018

8801SU.FG.023 +

801SU.FG.010

801SU.FG.014

801SU.FG.018

801SU.FG.023 +

6801SU.FG.016

6801SU.FG.018

6801SU.FG.023 +

5801SU.FG.016

5801SU.FG.018

5801SU.FG.023 +



8805SU.FG.014

805SU.FG.014



6813SU.FG.014



830SU.FG.009



8835SU.FG.010

8835SU.FG.012

835SU.FG.010

835SU.FG.012

6835SU.FG.010

6835SU.FG.012



5855SU.FGSS.016

8855SU.FG.016

8855SU.FG.018

8855SU.FG.025 ▮

6855SU.FG.016

6855SU.FG.018

6855SU.FG.025 ▮

5855SU.FG.016

5855SU.FG.018

5855SU.FG.025 ▮



5856SU.FGSS.016

5856SU.FGSS.018

8856SU.FG.012 +

8856SU.FG.014

8856SU.FG.016

8856SU.FG.018

8856SU.FG.021 +

856SU.FG.016

856SU.FG.018

856SU.FG.021 +

6856SU.FG.012 +

6856SU.FG.014

6856SU.FG.016

6856SU.FG.018

6856SU.FG.021 +

5856SU.FG.012 +

5856SU.FG.014

5856SU.FG.016

5856SU.FG.018

5856SU.FG.021 +



8856LSU.FG.016 +

8856LSU.FG.018 +

8856LSU.FG.020 +

856LSU.FG.016 +

6856LSU.FG.016 +

6856LSU.FG.018 +

6856LSU.FG.020 +

5856LSU.FG.016 +

5856LSU.FG.018 +



8850SU.FG.016 +

8850SU.FG.018 +

850SU.FG.012 +

850SU.FG.014 +

850SU.FG.016 +

6850SU.FG.012 +

6850SU.FG.014 +

6850SU.FG.016 +

6850SU.FG.018 +

5850SU.FG.012 +

5850SU.FG.014 +

5850SU.FG.016 +

5850SU.FG.018 +



8845KRSU.FG.016

8845KRSU.FG.018

8845KRSU.FG.025

6845KRSU.FG.016

6845KRSU.FG.018



8847KRSU.FG.016

8847KRSU.FG.018

847KRSU.FG.016

847KRSU.FG.018

6847KRSU.FG.016

6847KRSU.FG.018

5847KRSU.FG.016

5847KRSU.FG.018



8847SU.FG.014

8847SU.FG.016

8847SU.FG.018

6847SU.FG.014

6847SU.FG.016

6847SU.FG.018

5847SU.FG.014

5847SU.FG.016

5847SU.FG.018



8878KSU.FG.012 +

8878KSU.FG.014

8878KSU.FG.016

8878KSU.FG.018

878KSU.FG.012 +

878KSU.FG.014

6878KSU.FG.012 +

6878KSU.FG.014

6878KSU.FG.016

6878KSU.FG.018

5878KSU.FG.016

5878KSU.FG.018



8862SU.FG.010 +

8862SU.FG.012 +

8862SU.FG.014

862SU.FG.010 +

862SU.FG.012 +

862SU.FG.014

6862SU.FG.012 +

6862SU.FG.014



8888SU.FG.012 +

888SU.FG.012 +



8863SU.FG.012 +

863SU.FG.012 +

6863SU.FG.012 +



885SU.FG.012 +

6885SU.FG.012 +



5909SU.FG.040



8379SU.FG.018

8379SU.FG.023 +

379SU.FG.018

379SU.FG.023 +

6379SU.FG.018

6379SU.FG.023 +

5379SU.FG.018

5379SU.FG.023 +



8368SU.FG.016

8368SU.FG.023 +

368SU.FG.016

368SU.FG.023 +

6368SU.FG.016

6368SU.FG.023 +

5368SU.FG.016

5368SU.FG.023 +



8811SU.FG.033

8811SU.FG.037

6811SU.FG.033

6811SU.FG.037

5811SU.FG.033

5811SU.FG.037



8274SU.FG.016



8392SU.FG.016



KS0SU.FG.010 +



KS1SCSU.FG.012 +



KS2SCSU.FG.014



8858SU.FG.014 +

858SU.FG.014 +

6858SU.FG.014 +

5858SU.FG.014 +



8859SU.FG.010 +

8859SU.FG.016 +

859SU.FG.010 +

6859SU.FG.016 +

859SU.FGSL.010

Komet USA LLC
3042 Southcross Blvd, Suite 101
Rock Hill, SC 29730

Phone 888-566-3887
Fax 800-223-7485
info@kometusa.com
www.kometusa.com

single
patient use

Practice made perfect.

© Komet USA LLC · 02/2020 · 1001031\v001